Shannon Way, Tewkesbury, Gloucestershire. GL20 8ND United Kingdom Tel: +44 (0)1684 292 333 Fax: +44 (0)1684 297 929 187 Northpointe Blvd, Suite 105 Freeport, PA 16229 United States of America Tel: +1 724-472-4100 Fax: +1 724-472-4101 Tomson Centre 188 Zhang Yang Rd., B1602 Pudong New Area, Shanghai, Postal code: 200122 CHINA Tel/Fax: +86 21 587 97659 SCMC House 16/6 Vishal Nagar Pimpale Nilakh, Wakad, Pune PIN 411027 INDIA Tel: +91 827 506 5441





# APPLICATION NOTE

#### 1. Changes:

Version 1.0 Completed draft 23<sup>rd</sup> Jan 2019

#### 2. Introduction

 $\underline{FRAME} = 114$  enables the transformation for a 3 Arm Delta Robot. This allows you to define the end position of the wrist in X, Y, Z, and wrist angle (relative to the Z axis). The FRAME 114 transformation is available for MC4XX controllers (including MC464, EURO408 and MC405) with the Robotic Function FEC (feature enable code) fitted.

FRAME 114 is similar to FRAME 14 but with some additions. These additions include the inclusion of several optional parameters in the table data defining the robot. Further, while frame 14 is inflexible in that the robot parameters must be programmed in Micrometres, Frame 114 allows the robot to be defined in any desired unit, which will then be the units in which the robot is programmed.

#### 3. Axis order

The axes within the allocated FRAME\_GROUP have the following functions:

First Axis	Cartesian X axis
Second Axis	Cartesian Y axis
Third Axis	Cartesian Z axis
Fourth Axis	End effector rotation axis



#### 4. Diagram

The following pictures show the types of machine that can be controlled using this FRAME:



Once the FRAME is enabled set the UNITS to FRAME\_SCALE so that the Cartesian movements use the same scale as that used in the table data. So, if the TABLE data is programmed in mm then when UNITS is set to FRAME\_SCALE then the robot can be programmed in mm.

Set the UNITS on the rotational (wrist) axes to FRAME\_SCALE \* 2PI / 360 so that they are programmed in degrees. You can of course set UNITS for all axis to any suitable scale.

#### 5. Datuming

When datuming, (homing) it is necessary for each of the arms of the delta robot to be in the horizontal position (i.e. with the joints at either end of the link at the same level with respect to the robot's base).





## 6. Parameterization

13 parameters needed to be set into a sequence of TABLE values for use by the transformation mathematics. The base TABLE position is set using FRAME\_GROUP.

Table data	0	Top radius to joint (R1)
	1	Wrist radius to joint (R2)
	2	Upper arm length (L1)
	3	Lower arm length (L2)
	4	Edges per radian
	5	Angle of rotation in radians (Rotation)
	6	Linkx (optional with 4 or 5 axis)
	7	Linky (optional with 4 or 5 axis)
	8	Linkz (optional with 4 or 5 axis)
	9	Encoder edges/radian (optional Z rotation)
	10	Encoder edges/radian (optional Y rotation)

Note that R1, R2, L1, L2 should all be defined in the same UNITS that the robot will be programmed in. Also, note that the angle of rotation in radians, is the positive rotation of the X axis about the robot's base, with respect to the direction in which the axis 0 arm is pointing.

## 7. Example

BASE(0)

#### REPEAT

```
IF VR(1)=2 THEN
GOSUB frame_114 'Move the robot in Cartesian system
VR(1)=0
ELSEIF VR(1)=1 THEN
GOSUB frame_0 'Move the robot in Joint mode
VR(1)=0
ENDIF
UNTIL FALSE
STOP
```

frame\_0:

'configure for frame 0:



```
arm1_encoder = 16000 'Encoder counts per revolution
  arm1_gearbox = 360 '50
  arm1_edges_per_radian = arm1_encoder * arm1_gearbox / (2 * PI)
  arm1_edges_per_degree = arm1_encoder * arm1_gearbox / (360)
  arm2_encoder = 16000 'Encoder counts per revolution
  arm2_gearbox = 360 '50
  arm2_edges_per_radian = arm2_encoder * arm2_gearbox / (2 * Pl)
  arm2_edges_per_degree = arm2_encoder * arm2_gearbox / (360)
  arm3_encoder = 16000 'Encoder counts per revolution
  arm3_gearbox = 360 '50
  arm3_edges_per_radian = arm3_encoder * arm3_gearbox / (2 * Pl)
  arm3_edges_per_degree = arm3_encoder * arm3_gearbox / (360)
  rot_encoder = 200*32*16*80/360 'Encoder counts per revolution
  rot _gearbox = 360 '50
  rot _edges_per_radian = rot_encoder * rot _gearbox / (2 * PI)
  rot _edges_per_degree = rot _encoder * rot _gearbox / (360)
  FRAME = 0
  BASE(0)
  UNITS = arm1_edges_per_degree 'degrees
  AXIS_UNITS = UNITS
  RS_LIMIT=-1000
  FS_LIMIT=1000
  BASE(1)
  UNITS = arm2_edges_per_degree 'degrees
  AXIS_UNITS = UNITS
  RS_LIMIT=-1000
  FS LIMIT=1000
  BASE(2)
  UNITS = arm3_edges_per_mm 'mm
  AXIS_UNITS = UNITS
  RS_LIMIT=-1000
  FS_LIMIT=1000
  BASE(3)
  UNITS=rot_edges_per_degree
  AXIS_UNITS = UNITS
  RS LIMIT=-360
  FS_LIMIT=360
  FOR ax=0 TO 3
    SPEED AXIS(ax)=15
    ACC(200) AXIS(ax)
  NEXT ax
RETURN
```

frame\_114:



```
TABLE(0,400)
TABLE(0, 115)
TABLE(1, 44)
TABLE(2, 260)
TABLE(3, 715)
TABLE(4, 131072)
TABLE(5, pi/4)
```

'configure for frame 114:

```
FRAME_GROUP(0,0,0,1,2,3)'Delta robot
FRAME = 114
BASE(0)
UNITS = FRAME_SCALE 'mm
RS_LIMIT=-1000
FS_LIMIT=1000
BASE(1)
UNITS = FRAME_SCALE 'mm
RS_LIMIT=-1000
FS_LIMIT=1000
BASE(2)
UNITS = FRAME_SCALE 'mm
RS_LIMIT=-1000
FS_LIMIT=1000
BASE(3)
UNITS =( FRAME_SCALE * 2 * PI) / (360)'degrees
RS_LIMIT=-360
FS_LIMIT=360
FOR ax=0 TO 3
```

```
SPEED AXIS(ax)=10
ACC(200) AXIS(ax)
NEXT ax
```

#### RETURN